

# SERVER®

• E S S E N T I A L S •

## HoldCold™ Jars, Round DISCONTINUE NOTICE

PRODUCT UPDATE BULLETIN

### EFFECTIVE | JULY 5, 2022

- HoldCold™ Round Jars become obsolete
- Recommended cooling alternative; ice packs

### DETAILS

Due to supply chain disruption and low run sizes, HoldCold Jars have become obsolete; stock on-hand has run out. A chilled dispenser with alternate cooling method includes (2) ice packs for up to four hours of cooling between change-outs.

### PARTS AFFECTED

94055	HoldCold Jar, Round, Black
94056	HoldCold Jar, Round, White
94057	Kit of (1) Black & (1) White

### COMPLETE UNITS WITH AFFECTED PARTS

94111	CSM	Chilled Dispenser w/(1) 94055
94113	CSM	Chilled Dispenser w/(1) 94055 & 94056



### HOLDCOLD, ROUND JARS



### DISCONTINUED

CSM 94111 & 94113  
Cooled with HoldCold Jar(s)



### ALTERNATE COOLING METHOD

CSM 94000  
Cools with (2) 100515 Ice Packs  
and includes 94009 Stainless Steel Jar